

No: 10325-000

EIP 2: Irrigation Pump for Endoscopic P



- Good visibility of the target tissue
- EIP can be activated during intestinal irrigations either simultaneously or supplementary to diagnosis and therapy

Mains connection

Power frequency	50 / 60 Hz
Potential equalization connection	yes

Additional product information

Packaging unit	1 Piece
----------------	---------

Standards

Classification acc. to 93/42 EWG	Ila
Protection category acc. to EN 60 601-1	I
Type acc. to EN 60 601-1	CF